
# ecorun2.0

Teelicht-Stirling-Motor / Candle Stirling Engine

Before running the engine please read this instruction manual carefully!

Attention: The Stirling Engine ecorun2.0 is not a toy and not suitable for children.

### Safety Instructions

When running this engine use the same safety precautions as for open fire and flames. Never leave the engine running or the candle burning un-supervised. Never let a stopped engine heat for more than 3 minutes. Don't touch the engine's bottom cover - burn injury hazard.

Exergia is not responsible for any damages or injury occurring during the operation of this engine.

### Introduction

In the early days of the Industrial Revolution, the Edinburgh brothers James and Robert Stirling started working in 1815 on a safer alternative to the established steam engines of the day. By using air instead of steam as the working fluid, the design avoided the then common risk of boiler explosion. The first patent in 1816.

The Stirling engine is a vastly different heat engine to the internal-combustion (diesel or petrol) engine in your car. It has the potential to use any heat source as "fuel" including waste heat and solar energy. As the limits of fossil fuel resources became obvious, scientists and engineers revisited the Stirling principle as a way to exploit alternative energy sources. During the 1980s Prof. I. Kolin (University of Zagreb) and others developed engines based on an adapted Stirling principle capable of operating with a temperature differential of less than 50 degrees Celsius, thus proving the feasibility of exploiting low temperature heat sources to power modern Stirling engine designs.

The ecorun2.0 Stirling engine is an extraordinary engine that demonstrates these concepts with the conversion of medium temperature heat into mechanical work.

### Design

The ecorun2.0 Stirling engine consists of five main units:

#### X Main working space

'Can' comprising top and bottom housing covers connected by a white cardboard ring

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A graphite piston running in an aluminium power cylinder directly connected to the main working space.

**X Displacer piston**

A foam disc inside the main working space.

**X Crankshaft, con rods and propeller**

Two discs with pins connected via the main axle built the crankshaft. Two rods connect the crankshaft to the displacer and power piston – the phase difference between the sinusoidal movement of the pistons is about 90 degrees.

**X Rack with candle**

## Function

The bottom plate of the engine is heated up by the candle in the rack. When the displacer piston moves up air flows through it down to the bottom plate and is warmed up. On its way back to the top plate of the engine hot air streams through the displacer piston and is cooled by the aluminium top plate. The periodic change of the air's temperature causes a variation of its pressure, because it is encapsulated in the engine's housing. This results in a periodic force on the working piston on top of the engine. This force is transmitted by the connecting rod and crankshaft to the propeller (flywheel) and converts the up-and-down movement of the working piston into a rotational movement. As the expansion of the hot air through the power piston gives more work than it is spent by compression of "cold" air the surplus is the useful work on the working shaft. The displacer itself gets the power for its movement from the rotating propeller via a second connecting rod.

## Operating Instructions

Before starting, make sure that you're set up on a fireproof base in a dry and draft-free environment.

Put the tea light candle into the rack, light the candle and fit the engine centred on top of the rack over the burning candle. After a minute or two, flip the prop in a clockwise direction (assuming you are looking toward the prop). The ecorun2.0 engine should start immediately and will run for more than three hours - the burning time of the candle.

Please note that the size of the candle flame may vary during burning resulting in the engine's speed varying anywhere between 200 and 600 rpm.

The candle flame must never touch the surface of the

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lower housing cover. If necessary, shorten the candle's wick with a pair of scissors.

Never use other heat source than the delivered combination of rack with a food warmer candle. This engine is a demonstration model and not made for continuous use!

## Maintenance and Troubleshooting

Your engine has been tested and should run out of the box without problems. If you can't get your engine to operate properly or not run at all, go through the following check list step by step:

### X Engine heating

To ensure sufficient heating, the candle flame needs to be at least 10 to 15 millimetres high. If necessary, replace the candle.

A draft-free environment is essential to ensure optimal heat transfer from the candle flame to the bottom housing cover.

### X Lubrication

The engine has five bearing bushes: two bushes for the main axle in the support plates, two bushes in the connecting rods and one in the top cover to guide the displacer rod. From time to time apply a tiny drop of machine oil to each pin or axle where it goes through the bushes. Make sure that the oil goes into the bearing. **DO NOT** get any oil on the power piston or cylinder!

### X Free running piston

If any moisture has entered the engine's working space, water may have condensed in the gap between the power piston and cylinder, causing the piston to stick. Disassemble the power piston together with its con rod and dry it. Leave it disassembled for a while so any residual moisture may evaporate.

### X Connection con rod and piston - spare part

If there is any need for replacement of the small flexible tubes connecting the con rods and related pistons use the spare part included to the candle's foil bag. From this tube cut 7 mm lengths and carefully push this piece about 3 mm over the bare ends of each pin.

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