

zu beziehen bei  
sold by  
[www.conatex.com](http://www.conatex.com)



# DESCRIPTION

zu beziehen bei  
sold by  
[www.conatex.com](http://www.conatex.com)



CL-01 is a 150W cold light generator with a high efficiency 12 V halogen dichroic lamp. The brightness can be adjusted using a remote control panel.

## SPECIFICATIONS

Halogen lamp:	150W, 12 V, 50 Hz, socket 6,35 (OSRAM HLX 64627 or equivalent)
Fuse:	2 A T
Power supply:	230 V

## SAFETY PRECAUTIONS AND MAINTENANCE

### Changing the halogen lamp

Before opening the illuminator, turn off the tension, unplug the power supply and wait for ten minutes until the lamp has cooled down. Unscrew the screw on the top of the illuminator and open the lid. Remove the old lamp by pulling it upwards and by carefully pulling the pins out of the socket. Replace with a new lamp and close the lid with the screw.

**Warning:** during use, the lamp becomes very hot. This is normal, and does not cause malfunction. Beware to keep flammable objects at, at least 0.5 m away from the light source. Do not obstruct the fan and the venting holes.

### Cleaning the instrument

Clean the generator with a soft dry cloth. Do not use solvents.

# RECOVERY AND RECYCLING

zu beziehen bei  
sold by  
[www.conatex.com](http://www.conatex.com)



Art. 13 Dlsg 25 July 2005 N°151. "According to directives 2002/95/EC, 2002/96/EC and 2003/108/EC relating to the reduction in the use of hazardous substances in electrical and electronic equipment and waste disposal."



The basket symbol on equipment or on its box indicates that the product at the end of its useful life should be collected separately from other waste.

The separate collection of this equipment at the end of its lifetime is organized and managed by the producer. The user will have to contact the manufacturer and follow the rules that he adopted for end-of-life equipment collection. The collection of the equipment for recycling, treatment and environmentally compatible disposal, helps to prevent possible adverse effects on the environment and health and promotes reuse and/or recycling of materials of the equipment. Improper disposal of the product involves the application of administrative penalties as provided by the laws in force.