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Fully automatic **GFL Water Stills models 2002, 2004, 2008 and 2012 with storage tank** produce highly-pure, bacteria and pyrogen free distillate with a very low conductivity (approx 2.3 $\mu\text{S}/\text{cm}$ at 25°C). The distillate is in conformity with DAB regulations and the regulations of many international pharmacopeia.

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- Main switch
- Pilot lamp Clean
- Operation pilot lamp
- Distillate withdrawal

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Before installation, please check whether contents of package are in good order and complete.

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Should you note any damages or have any reasons for complaint, please contact your supplier

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Contents –

1	Using the Water Still	Page	12
2	Warranty	Page	12
3	Before Initiation	Page	13
4	Transport, Location and Set-up of the Water Still	Page	13
5	Operation Voltage	Page	14
6	Water Connections	Page	14
7	Initiation	Page	16
8	Functional Description	Page	16
9	Maintenance, Support and Trouble Shooting	Page	17
10	Disposal of Used Units	Page	18
11	Technical Data	Page	19
12	Circuit Diagram	Page	20
13	Connection to the Mains	Page	24
14	Spare Parts List	Page	29
15	EC Declaration of Conformity	Page	32
16	Accessories	Page	33

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1. Using the Water Still

1.1 Intended Use

- >  GFL Water Stills have been tested by TÜV SÜD for compliance to the basic requirements of the laboratory equipment regulations EN 61010-1:2010 and EN 61010-2-010:2014 as well as to the general engineering standards. The manufacturing site is regularly controlled by the testing institution. The units may be marked with the depicted quality symbol of TÜV-SÜD Product Service GmbH.
- > In Mono Water Stills 2002 to 2012, electric heating elements serve to boil water and evaporate it. The steam is led through steam tubes into the storage tank where it condenses on a water-cooled cooling coil. The produced distillate drips into the storage tank. The storage tank level is electronically monitored, the Water Still's heating element(s) and cooling water supply will be switched off when the tank is full. Depending on the model, the distillate quantity per hour ranges between 2 litres (2002), 4 litres (2004), 8 litres (2008) and 12 litres (2012).
- > Depending on the tap water quality, the produced distillate has a conductivity of approx. 2.3 µS/cm at 25°C. To feed the Water Still, preferably use tap water of drinking quality. Please also observe the information in chapters 16.2 to 16.4 of these instructions on the possible use of filters and cartridges to pretreat the tap water.
- > The Water Stills must be operated within the user's field of vision.
- >  The information contained in these operating instructions must **by all means** be read and observed. Only then can a perfect operation of the **Water Still** be guaranteed. The units may only be installed and operated by persons who have made themselves familiar with these operating instructions.
- >  The accessible inner parts of the housing that can be reached after lifting the outer lid get strongly heated up during operation. These parts may only be touched after they have cooled down, or when wearing suitable safety gloves.

1.2 Improper Use

- > GFL Water Stills, operated in a laboratory, are **no Medical Devices**. They fall neither under national nor international Medical Device Directives and have to be used and applied accordingly.
- > The Water Still may not be used in potentially explosive surroundings.
- > The Water Still may neither be set up nor operated in laboratory areas with aggressive or corrosive ambient conditions.

2. Warranty

The warranty period for all laboratory apparatus and their accessories by GFL Gesellschaft für Labortechnik mbH is 24 months from the date of the invoice made out by GFL. The warranty period for repairs and modifications, carried out by GFL, as well as spare parts is 6 months from the date of the invoice made out by GFL. In order to identify defective units, we require both model and serial number on the nameplate at the back, left-hand side of the Water Still and, if applicable, a copy of the invoice.

A free-of charge guarantee repair will not be granted for defects that are due to improper installation and handling.

3. Before Initiation

Important information are marked in bold letters in these instructions, safety indications are additionally marked by the following warning symbols and mandatory signs.

Warning of hot liquids and vapour

Warning of hot surfaces


Read and observe the operating instructions

Warning of dangerous electrical voltage

General warning

Before maintenance and repair disconnect the unit all-pole from the electrical mains (pull the plug from the socket).

All free openings of the water connections on the right-hand side of the unit are taped. These transport protections. These protections have to be removed before initiation.

4. Transport, Location and Set-up of the Water Still

Protect yourself and the unit during transport and setup by working carefully and avoid danger of e.g. shifting or tilting the unit as well as risk of injury by lifting heavy loads.

Caution, Water Stills 2004 (net weight 21.8 kg), 2008 (net weight 35.8 kg) and 2012 (net weight 40.9 kg) must be lifted, carried and transported to the location by at least two persons.

The Water Still can be held between the four stands of the unit and lifted for setup.

The unit is suitable for both bench and wall mounting. It is only suitable for indoor use.

Table-top setting-up on solid, even and level surfaces only. Make sure to place the unit only on a watertight, temperature-resistant, non-flammable surface. The location must provide sufficient space as well as the necessary carrying capacity for the total weight of the unit (unit weight as per technical data, chapter 11 of this manual, plus weight of the filling).

For **setting-up on the wall** check the carrying capacity of the wall in connection with the total weight of the unit (appliance plus water filling, see Technical Data). Make sure to check the on-site fixing elements in order to guarantee safe hold of the Distiller. The fixing material must be suitable for the substrate. Use only tested and proven fixing material.

There are two keyholes for wall mounting at the back of the **Water Still**. Place two fixing screws (not included in the scope of supply) in the wall in the distance of the two keyholes and nest the still on the keyholes. The unit is delivered without fixing material.

Model / Distance of screws

2002 / 40 cm 2004 / 48 cm

2008 / 64 cm 2012 / 64 cm

The unit is not suitable for use in potentially explosive surroundings, e.g. during anaesthesia with inflammable gas or steam types!

5. Operation Voltage

The Water Still must be connected to the mains supply either by a correctly installed shock-proof socket, or, with a fixed connection, to an on-site main switch. **The Water Still is a protection class I electrical appliance, a connection to the earth conductor (PE) must be ensured.** For information on the required mains fuse please view Technical Data, chapter 11 of this manual.

The electrical connection must ensure an all-pole separation from the mains at any time. Shock-proof sockets or main switches to separate the Water Still from the mains must be installed in such a way as to ensure they can be clearly identified and are easily accessible at any time.

The mains connection cable must not touch any hot surfaces of the unit. It may not lead underneath the unit.

The main switch (I) of the Water Still must be **off** (position O). **The voltage on the nameplate (at the left-hand side of the unit) must be identical to the mains voltage. If they are identical, connect the unit to the mains.**

Please also refer to chapter 13 of these operating instructions "Connection to the Mains".

6. Water Connections

All water connections of the Water Still, except the distillate withdrawal tap on the front, are situated on the right-hand side of the unit.

Hoses for water inlet and outlet are not included in the standard scope of supply.

6.1 Distillate withdrawal ①

Distillate is withdrawn from the unit through the black plastic tap on the front of the unit. It can be opened in continuous or interval positions. A laboratory hose with an inner diameter of approx. 15 mm can be connected to the tap, the hose has to be secured from slipping off with a hose clip.

Caution, distilled water leaves the Water Still with a temperature of more than 50°C.

Danger of scaldings!

6.2 Drain boiler ②

Connect a ½" hose to drain the boiler when cleaning the unit and/or carrying out maintenance works.

Caution, the water leaves the boiler of the Water Still with a temperature of up to 100°C.

Danger of scaldings!

6.3 Cooling water outlet ③

Connect a $\frac{3}{4}$ " **temperature-resistant hose** to the cooling water outlet. It may have a maximum length of approx. 1.5 m. It must be led into an open drain on a lower level and must have a slope on its complete length.

Make sure the cooling water can drain without back draughts.


**Caution, cooling water leaves the Water Still
with a temperature of up to 70°C.
Danger of scaldings!**

6.4 Tap water inlet ④

The tap water inlet supplies water to the water still through a solenoid valve. The hose connection of the valve must be connected with the water mains through a $\frac{1}{2}$ " **pressure hose** and a blockable water connection (⑤, hand stop valve).

It is absolutely essential to use hose clips to secure both hose connections!


7. Initiation

7.1 Before initial starting

Before initial starting, the Water Still has to be **hand-filled with water**. To do so, remove the outer and left-hand inner lids (marked **L**) as well as the deflector beneath the inner lid. Caution, when working on the deflector, make sure to observe the correct position of the wire electrode in the container. It may neither be bent nor come into contact with the housing. The boiler must now be filled with water until the heating element(s) at the bottom of the container are below water level. The deflector as well as the two individually marked inner lids must be re-assembled and safely nested on the containers before initiation of the unit.

Marking L
front of the left-hand
inner lid (boiler)

Marking R
front of the right-hand
inner lid (storage tank)


7.2 Initiation

To initiate the Water Still after assembly of all power and water connections, open the hand-stop valve of the water mains supply (5) and switch on the main switch (**I**). **The first few litres of produced distillate should not be used. Discard the first two to three complete fillings of distilled water in the storage tank.**

To switch off the unit for longer down-times first switch the main switch (**I**) to position O, and disconnect the unit from the mains. Close the on-site hand-stop valves of the water supply, then empty and dry storage tank and boiler in order to avoid germ contamination.

8. Functional Description

I

III

GFL Water Stills 2002, 2004, 2008 and 2012 work automatically.

After switching on the Water Still's main switch the green pilot lamp in the main switch (**I**) as well as the right-hand, yellow operation pilot lamp (**III**) will light up. The built-in solenoid valve opens. Water flows through the cooling coil in the storage tank and a mechanical water level regulator that controls the water level in the boiler. Excess water that is not used to evaporate, runs off through the cooling water outlet. The heating elements are powered on

and bring the water in the alembic to the boil. A thermostatic **low water cut-off** protects the heating elements from dry running. The produced steam is led through a steam tube to the cooling coil where it condenses and drips as distillate into the storage tank.


The water level in the storage tank (right-hand container) is controlled through a wire electrode. Once the container is full, the Water Still will be switched off by an electronic controller. The solenoid valve interrupts the cooling water supply, the heating elements are switched off and the yellow pilot lamp (**III**) goes out.

Distillate is withdrawn through the black plastic tap (IV) at the front of the unit.
After withdrawal of distillate the unit will be switched on again automatically so that the storage tank is refilled.
Carbon dioxide is degassed through a vent on top of the unit.

9. Maintenance, Service and Trouble Shooting

**Caution! Prior to maintenance and service on the unit let the Water Still cool down!
Danger of burns!**

**Caution, before opening the Water Still as well as prior to cleaning separate the unit from the mains supply!
(Pull the plug or switch off main switch)
Danger of electrical shock!**

9.1 Descaling

Depending on the degree of hardness of the tap water, the Water Still must be cleaned regularly of scale formations. First remove the outer and the left-hand inner lids (L), then unscrew the deflector and remove it, too.

A suitable descaling agent is a mixture of 10% formic acid, 10% acetic acid and 80% distilled water. Fill the solvent into the boiler until it covers the topmost scale deposits and **heat up to maximum approx. 70°C**. Make absolutely certain not to bring the mixture to a boil. In order to do so, close the hand-stop valve of the water mains supply (5) and switch on the Water Still until the temperature is reached. After approx. 30 minutes drain the solvent and scale mixture through the drain cock "drain boiler" (2) and **rinse the boiler thoroughly with water several times**.

Commercial descaling agents, suitable for use with stainless steel, can also be used according to the manufacturer's instructions (e.g. rea-calc® of M/s CHEMOTEC GmbH, 63486 Bruchköbel, Germany).


**Never use any products containing hydrochloric acid!
These will lead to damages to heating element, boiler,
temperature sensor and the ducts of the screw connections.**

Restart the unit as described in chapter 7 **Initiation**.

After descalings, the first few litres of distillate should not be used as it might contain traces of evaporated descaling agent. Discard the first two to three complete fillings of distilled water in the storage tank.

9.2 Pilot Lamp Clean

9.2.1 Foam formation caused by polluted water in the boiler

I

II

Depending on the degree of impurities in the feed water as well as on the growing contamination of the water in the boiler, caused by the distillation process, the boiling water will foam up. Once the foam gets in contact with the electrode in the boiler (L), an electronic impurity detector will switch off the unit, and the **red pilot lamp Clean (II) will light up**. The boiler must now be emptied through the drain cock "drain boiler" (2) and be rinsed several times with clean water.

The operational interruption "Clean" is reset by switching the unit's main switch off and then on again.

The subsequent re-initiation is to be carried out as described in **chapter 7**.

The "Cleaning" function's purpose is to exchange polluted water in the boiler against clean water.

9.2.2 Back water in the boiler

If the water level in the boiler rises up to the wire electrode in the left-hand container (**L**), the electronic impurity detector (A1 in the circuit diagram) switches the unit off and **the red pilot lamp "Clean" (II)** lights up. This means that more water entered the unit than could flow off. Possible causes are tubes outside the unit that are not installed as described (see chapter 6), or scale formations that cause a blockage in the water outlet tube inside the unit. A defective quantity regulator in the outlet of the solenoid valve (see cooling water requirement, chapter 11) might also cause more water to flow into the unit than can flow off. For trouble shooting purposes, an additional function drawing of the internal water flow is available on request.


9.3 Re-initiation after low water

A low water cut-off (thermostatic over-temperature cut-out) protects the heating elements against running dry. In case of low water the Water Still is switched off by the temperature limiter.

To re-initiate, let **the unit cool down** and fill in water to cover the heating elements as described in chapter 7 **Initiation**.

The triggered low water cut-off must be reset. To do so, loosen the black cap nut at the bottom left-hand side of the Water Still. Inside the thread a small white plastic pin can be seen that has to be gently pressed inside (e.g. with a pen) until a clicking sound can be heard.

Only then switch the unit back on again.

The **GFL Water Still** is made of first-class material. Nevertheless, the unit should only be subjected to mechanical strain within sensible limits.

Make sure that no fluids come into contact with cable connections or the electrical parts inside the unit!

The powder-coated surfaces of the unit can be cleaned with mild detergents, if necessary.

Servicings, repairs or modifications must be carried out according to the General Engineering Standards by a competent electrician.

Only original spare parts must be used. Always demand a detailed confirmation by the person in charge (company, date, signature) on the kind and volume of the tasks carried out.

10. Disposal of Old Units

GFL will take responsibility, within the scope of the legal directives, for an environmentally sound handling and disposal of all used GFL units as of the production year 1995 that are returned to us **free of charge** and will have it materially recycled. Before the unit is returned, a legally binding declaration must be provided from the sender confirming that the unit is free from harmful and/or hazardous contaminations as well as from hazardous substances caused by the previous use of the unit.

GFL laboratory apparatus are exclusively designed for industrial use and may not be disposed of through public waste disposal authorities.

EAR Registration Number WEEE-ID.NO.DE 67770231


11. Technical Data

Exterior dimensions (W x D x H)	Model 2002	530 x 280 x 455 mm	Model 2008	780 x 405 x 575 mm
	Model 2004	615 x 320 x 495 mm	Model 2012	780 x 405 x 705 mm
Storage tank	Model 2002	4 litres distillate	Model 2008	16 litres distillate
	Model 2004	8 litres distillate	Model 2012	24 litres distillate
Distillation capacity	Model 2002	2 l / h distillate	Model 2008	8 l / h distillate
	Model 2004	4 l / h distillate	Model 2012	12 l / h distillate
Distillate quality	approx 2.3 µS/cm at 25°C in conformity with DAB, bacteria and pyrogen free, low gas content.		The conductivity of the distilled water is directly related to the chemical composition of the raw water. Components of the raw water having the same or a lower evaporation point than water may deteriorate the conductivity.	
Cooling water requirement	Model 2002	30 litres / hour	Model 2008	72 litres / hour
	Model 2004	48 litres / hour	Model 2012	198 litres / hour
Low water cut-off	Electromechanical temperature limiter, with capillary tube sensor. Switch-off temperature 135°C -15 K			
Water pressure min. / max.	> 3 bar / 7 bar > 43.5 psi / 101.5 psi			
Electrical connection Mains supply / On-site fuse / Internal mains fuse	Model 2002	230V +/-10%, 50/60 Hz, 1.5 kW shock-proof plug on-site 16 A / internal 10A T	Model 2008	400V / 3 / N / PE +/-10% 50/60 Hz, 6.0 kW Mains connection cable for permanent connection to the mains on-site 3x16 A
	Model 2004	230V +/-10%, 50/60 Hz, 3.0 kW shock-proof plug on-site 16 A / internal 16A T	Model 2012	220V / 3 / PE +/-10% 50/60 Hz, 9.0 kW Mains connection cable for permanent connection to the mains on-site 3x25 A /
Caution! Mains voltage deviations, even with the mentioned limits, influence the quantity of distillate produced.	Model 2008	230V +/-10% 50/60 Hz, 6.0 kW Connection box for permanent connection to the mains on-site 35 A /	Model 2012	400V / 3 / N / PE +/-10% 50/60 Hz, 9.0 kW Mains connection cable for permanent connection to the mains on-site 3x16 A
	Model 2008	220V / 3 / PE +/-10% 50/60 Hz, 6.0 kW Mains connection cable for permanent connection to the mains on-site 3 x 16 A		
Protection class / - type	I / IP20			
Surrounding conditions	Height above MSL			
	up to 2000 m MSL			
Ambient temperature	+10 °C to +40 °C			
Humidity	max. 80% rel. humidity to 31°C, decreasing to 40% rel. humidity at 40°C			
Weight	Model 2002	16 kg / 22 kg	Model 2008	35.8 kg / 56.5 kg
	Model 2004	21.8 kg / 32.4 kg	Model 2012	40.9 kg / 70.4 kg
Net / with water filling				

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
s. Plaque d'identification

13. Connection to the Mains

Examples for Connection to the mains supply pages 25 to 28

The electrical connection must ensure an all-pole separation from the mains supply. Installed assembly parts, such as shock-proof sockets or main switches, must be installed so as to ensure clear identification and they must be within easy reach in an emergency case. **Water Stills 2008 and 2012 must have a permanent connection with the mains. This may only be implemented through an on-site main switch or through a CEE plug in conformity with IEC standard 60309-2 (see examples for connection to the mains in chapter 13.2).**

Colour decoding of the individual leads of the mains connection cables for models 2008 and 2012

Colour decoding	Mains supply		Mains supply	
	400V / 3 / N / PE	50/60 Hz	220V / 3 / PE	50/60 Hz
ge/gr – gelb/grün – yellow/green	PE		PE	
bl – blau – blue	N			
sw – schwarz – black	L1		L1	
br – braun – brown	L2		L2	
gra – grau – grey	L3		L3	

Water Stills can be supplied in different versions for connection to different mains supplies.

Models 2002 and 2004 for connection to 230V (see nameplate) can be connected **to all mains supplies with 220V or 230V.**

Models 2008 and 2012 for connection to 220/3 (see nameplate) can **only** be connected to mains supplies with **220V / 3 / PE.**

Models 2008 and 2012 for connection to 400/3 (see nameplate) can **only** be connected to mains supplies with **400V / 3 / N / PE.**

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Model	Power	Power consumption at mains voltage*	Mains fuse (F2–F4)
2002	1,5 kW	6.5 Amp at 230 V	10 Amp
2004	3,0 kW	13.0 Amp at 230 V	16 Amp
2008	6,0 kW	26.1 Amp at 230 V	35 Amp
	6,0 kW	15.8 Amp at 220V / 3 / PE	16 Amp
	6,0 kW	8.7 Amp at 400V / 3 / N / PE	10 Amp
2012	9,0 kW	23.6 Amp at 220V / 3 / PE	25 Amp
	9,0 kW	13,0 A bei 400V / 3 / N / PE	16 Amp

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

Mains connection through shock-proof plug CEE 7/7 (pre-assembled), all-pole separable from the mains.

13.1.2 Type 2002 et 2004 en 230V sur une alimentation 400V / 3 / N / PE / 50/60 Hz Models 2002 and 2004 for 230V for a mains supply of 400V / 3 / N / PE / 50/60 Hz (s. Plaque d'identification/ see nameplate)

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

Mains connection through on-site switch S4, all-pole separable from the mains. Connection cable not in scope of supply.

13.1.4 Type 2008 en 230/1 à une alimentation 400V / 3 / N / PE / 50/60 Hz Models 2008 for 230V for a mains supply of 400V / 3 / N / PE / 50/60 Hz (s. Plaque d'identification / see nameplate)

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
Mains connection through CEE plug in conformity with IEC standard 60309-2, all-pole separable from the mains. Connection cable and CEE plug not in scope of supply.

13.1.5 Type 2008 et 2012 en réseau électrique 220/3 220V / 3 / PE / 50/60 Hz
Models 2008 and 2012 for 220/3 for a mains supply of 220V / 3 / PE / 50/60 Hz
 (s. Plaque d'identification / see nameplate)

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
Mains connection through on-site switch S4, all-pole separable from the mains.

13.1.6 Type 2008 et 2012 en alimentation 220/3 220V / 3 / PE / 50/60 Hz
Models 2008 and 2012 for 220/3 for a mains supply of 220V / 3 / PE / 50/60 Hz
 (s. Plaque d'identification / see nameplate)

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Mains connection through CEE plug in conformity with IEC standard 60309-2 (plug not pre-assembled), all-pole separable from the mains.

13.1.7 Type 2008 et 2012 en 400/3 sur une alimentation 400V / 3 / N / PE 50 / 60Hz
Models 2008 and 2012 for 400/3 (see nameplate)
for a mains supply of 400V / 3 / N / PE 50/60 Hz
 (s. Plaque d'identification / see nameplate)

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


declare that the below stated **Water Stills** models:

2002, 2004, 2008 and 2012

with the technical data:

**230V, 50 / 60 Hz
1.5 kW (2002)
3.0 kW (2004)**

**220V / 3 / PE, 50 / 60 Hz or
400V / 3 / N / PE, 50 / 60 Hz
6.0 kW (2008)
9.0 kW (2012)**

are in conformity with the following EC Directives:

I	2014/35/EU	(Low Voltage Directive)
II	2014/30/EU	(EMC Directive)
III	2011/65/EU	(RoHS 2 Directive)

For conformity **with I** the following standard was applied:

**EN 61010-1:2010
EN 61010-2-010:2014**

For conformity **with II** the following standard was applied:

EN 61326-1:2013

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Andreas Degmayr
Managing director

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16.1 Separate Water Supply,

to feed the boiler with softened or desalinated water (pressure > 1 bar / 14.5 psi) and the cooling coil with phosphatised or normal tap water (pressure > 3 bar / 43.5 psi). **Accessory must be installed in our works.** When a separate water supply is installed, the hourly capacity of distilled water is reduced by approx. 10-15 %.

Separate Water Supply

Order No 2901

Without illustration

Factory installation into models **2002 to 2012 possible.** When the storage is full, the inlet of pretreated water is not automatically switched off.

Separate Water Supply with
solenoid valve

Order No 2909

Without illustration

Factory installation into models **2004 to 2012 possible.** When the storage tank is full, the solenoid valve switches off pretreated water automatically.

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Dechlorite Filter,

eliminates chlorine particles in tap water added by the local waterworks. Complete with connections for pressure hose ½ inch, with first filling.

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Phosphate cartridge,

prevents scale formation in the condenser by phosphatising of tap water. Can be used from 4 to 15°dH (German hardness), equivalent to approx. 0.7 to 2.7 mMol/l. With connection for pressure hose ½ inch, with first filling.

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Pre-Filter 1 µm,

for pre-cleaning the tap water, and to protect the unit from premature contamination. Complete with connections for pressure hose ½" (inner diameter 12.7 mm), including filter candle. The candle should be replaced at least every six months.

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Hose Set,

consisting of hoses for water inlet and outlet (length 1.5 m), including hose clips.

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
Wall bracket,

for one filter or combinations of two or three filters, including sleeves to connect the filters and screws to fix the filters to the wall bracket.

A data sheet on possible connection variants of articles 2904, 2906 and 2912 is available on request.

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Wall bracket
for one filter
Order no 2912

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Wall bracket
for two filters
Order no 2922

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Wall bracket
for three filters
Order no 2923

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Level Control Switch

When connecting an external storage tank (not included in the scope of supply) to the Water Still, it is recommended to equip the unit with a Level Control Switch 2910 which controls the water level in the external storage tank and switches off power and water when the external storage tank is full. When using an external storage tank in connection with Level Control Switch 2910, the internal storage tank is only usable to a limited extent due to different levels of the internal and external tanks. This accessory must be installed in our works.

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Level Control Switch
Order No 2910