WPA S1200⁺

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EXTRAIT DE LA NOTICE ORIGINALE

WPA S1200+ 80-3007-20

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
Key	Action
On/Standby Key 	Turns the instrument on or back to standby.
Arrow Keys 	Use the four arrow keys to navigate around the display and select the required setting for the active (highlighted) option.
View Options 	Used to select items in a menu or to enter parameters and text descriptions. Repeated key presses can be used to cycle through lower case, upper case, and numbers. Wait 1 second before entering the next character. Use the "C" clear key to backspace and the 1 key to enter a space
Alphanumeric Keys 	Use to enter parameters and to write text descriptions where appropriate or required. Use repeated key presses to cycle through lower case, upper case, and numbers. Wait 1 second before entering the next character. Use the "Esc" key to backspace and to enter a space
Escape/Cancel 	Escape from a selection and return to the previous folder. Stop making measurements.
OK 	Confirm parameters/selection.
Set Reference 	Set reference to 0.000A or 100%T on a reference solution at the current wavelength in the measurement mode selected. When in scan mode, do a reference scan.
Take a Measurement 	Take a measurement.


Software Style

The user interface is built around having folders of files which are displayed on the home page when the instrument is switched on. Different folders are numbered and opened by using the associated number key on the keypad.

Folder	Keypad Number	Description
 Applications	1	Single wavelength, Concentration, Wavelength scan, Kinetics, Standard Curve, Multiple wavelengths and Ratio.
 Favourites	2	Saved User selected and configured methods (maximum of 9 methods).
 Methods	3	Sub folder selection for User selected and configured methods (9 sub folders each, containing a maximum of 9 methods).
 Utilities	4	Instrument set up (date, time, language, data output, preferences and display settings).
 Life Science	5	Proteins and Cell counting.

The Applications Folder

Function	Keypad Number	Description
 Single Wavelength	1	Absorbance or %T (transmission) at a single user defined wavelength.
 Concentration	2	Concentration measurement at a single wavelength based on a simple Factor entered or calculated from a single standard.
 Wavescan	3	Wavelength scan between two user-defined wavelengths. Range 330-950 nm, with user configurable peak finding function.
 Kinetics	4	Absorbance versus time measurements either rate or end value based.
 Standard Curve	5	Generation of calibration curve by measuring standards at a single wavelength.
 Multi Wavelength	6	Absorbance or %T (transmission) at up to 5 user defined wavelengths.
 Absorbance Ratio	7	Ratio of absorbance values at two user specified wavelengths.

Options

Within each application the user has the possibility to select various options that define the way results are treated. If not using a stored method, it is advisable to check that these Options have been appropriately set for your experiment when coming to the instrument. Note that setting the “History” parameter to on (see Preferences later) will cause the instrument


to store its last settings. If the “History” parameter is turned off, all parameters and options will return to their default settings when you leave that application. (Unless it has been saved as a method).

Applications

Single wavelength

This makes simple absorbance (A) and % transmission (%T) measurements on samples, measuring the amount of light that has passed through a sample relative to a reference (this can be air).

The procedure is as follows:

Step 1

Set wavelength by using keypad numbers or left and right arrows. Press the down arrow key.

Step 2

Select the mode, *Absorbance* or *%T*, using the left and right arrows. Press the down arrow key.

Step 3


Select whether or not to adjust the name of your sample before each new run, using left and arrow keys to select Yes or No.

Step 4

To enter the results screen with the selected

parameters press OK  or Cancel the selections and return to the Applications Folder by pressing

Cancel .


Step 5

Insert the reference. Press the reference key . This will be used for all subsequent samples until changed.

Step 6

Insert sample and press .

Repeat step 6 for all samples.

Results Screen

The result at the selected wavelength is displayed on screen.

Use the left and right arrows to move the cursor and display the value at the cursor position (+/- 15nm from set wavelength).

Press Cancel  to return to the Applications Folder.

Press  to display available options which are described below.

- ① Parameters...
 - ② Print
 - ③ Abs/%T
 - ④ Print Graph
-
- ⑦ Sample Number...
 - ⑧ Save Method...
 - ⑨ Auto-Print

Options (select using key pad numbers)

1. Return to parameters screen (Step 1 above).
2. Print result via selected method.
3. Toggle between Absorbance and %T mode.
4. Print graph – greyed out if no data are available.
7. Sample number – add a prefix to the sample number and reset the incrementing number to the desired value.
8. Save method – use the left and right arrows to select a folder to store in (Favourites/Methods 1-9), press the down arrow and enter name.
9. Auto-print – toggles auto-print on/off.

Exit options by pressing  or wait.